# Geistlich Bio-Gide® Shape

ENGLISH

## Geistlich Bio-Gide® Shape

Pre-shaped, Resorbable Bilayer Collagen Membrane for Guided Tissue and Bone Regeneration

Geistlich Bio-Gide® Shape is a pre-shaped biocompatible, sterile resorbable bilayer collagen membrane for guided tissue and bone regeneration.

### COMPOSITION

Geistlich Bio-Gide® Shape is a pure collagen membrane obtained by standardized, controlled manufacturing processes. The collagen is extracted from veterinary certified pigs and is carefully purified to avoid antigenic reactions. Geistlich Bio-Gide® Shape is sterilized in double blisters by  $\gamma$ -irradiation. Geistlich Bio-Gide® Shape has a bilayer structure: The porous surface – facing the bone – will allow the ingrowth of the bone forming cells. The smooth surface – facing the soft tissue – will prevent the ingrowth of fibrous tissue into the bone defect. The membrane is made of collagen without further cross-linking.

### PROPERTIES / ACTION

The low antigenicity and excellent biocompatibility support the use of Geistlich Bio-Gide® Shape in dental surgery. Animal studies have shown that the collagen is resorbed in bone cavities within 24 weeks. The fibrous microstructure of Geistlich Bio-Gide® Shape is hydrophilic. The membrane retains its structural integrity even when wet. Fixation by sutures or pins is possible. Geistlich Bio-Gide® Shape has a smooth, dense surface which facilitates cutting the membrane. The smooth, dense surface also prolongs the water uptake, application is facilitated. The coherent collagen fibres swell and form a unified basic tissue structure. As a result, the adaption to the bone wall and the complete closure of the bony defects is easily achieved. Inflammatory reactions have not been observed, but due to the special physical properties and elongated resorption time a certain inflammatory reaction after placing the membrane cannot be excluded.

### INDICATIONS FOR USE

Geistlich Bio-Gide® Shape is indicated for:

- augmentation around implants placed in immediate extraction sockets.
- augmentation around implants placed in delayed extraction sockets.
- localized ridge augmentation for later implantation.
- alveolar ridge reconstruction for prosthetic treatment.
- filling of bone defects after root resection, cystectomy, removal of retained teeth.
- guided bone regeneration in dehiscence defects.

### CONTRAINDICATIONS

Geistlich Bio-Gide® Shape should not be placed where active infection exists. Before placement, the surgeon should be confident that any active or recent infection has been properly treated.

Geistlich Bio-Gide® Shape should not be used in patients with a known allergy to collagen.

### DIRECTIONS FOR USE

The general principles of sterile handling and patient medication must be followed when using Geistlich Bio-Gide® Shape:

- The tooth is extracted and basic surgical procedures (e.g., curettage, de-epithelization) are undertaken.
- The long strap of the Geistlich Bio-Gide® Shape should be placed to cover the bony defect of the socket, if present. Thereby the smooth surface should face towards the soft tissue, and the rough side should face the defect.
- Place the membrane without pre-wetting within the extraction socket and hold in place with moderate pressure. Adhesion to the bone surface is facilitated by the hydrophilic properties of the membrane.
- The bone defect is filled with a space-making material, e.g., autogenous bone, Geistlich Bio-Oss®, Geistlich Bio-Oss Collagen®, etc. The defect should not be overfilled.
- Complete penetration of the membrane by blood and exudate allows excellent adaptation and adhesion of the membrane to the defect.
- To close the extraction socket, both lateral wings of the Geistlich Bio-Gide® Shape membrane should be tucked underneath the mucosa, located mesially and distally of the socket. The upper wing should be tucked lingually / palatally in-between the soft tissue and bone wall.
- Due to the high tensile strength of the membrane, fixation is possible. Fixation of the membrane (e.g. mattress suture) on top may be indicated to avoid its displacement due to shear loading or mobilization.
- Adhesion to the bone surface is facilitated by the hydrophilic properties of the membrane.
- Complete penetration of the membrane by blood and exudate allows excellent adaptation and adhesion of the membrane to bony structures.
- Salivary and other contamination to the material and surgical site should be minimized to avoid bacterial contamination.
- In case of membrane exposure (e.g. soft tissue dehiscence) during the healing phase the resorption time may be accelerated.

### POST-OPERATIVE CARE

- The patient should be monitored closely.
- If the membrane becomes exposed, the dehiscence usually heals by itself within several weeks. Membrane removal is usually not necessary. However, to minimize bacterial contamination rinsing with bactericidal solutions is recommended.
- In the event that early membrane removal is necessary, the tissues adjacent to the membrane should be anesthetized with a local anesthetic. An incision should then be made immediately adjacent to the residual membrane. Following careful reflection of the surrounding tissue, the remaining portion of the membrane can be excised and the area curetted to remove any inflamed or infected tissue.
- To allow for undisturbed bone regeneration underneath the membrane surgical reentry should not be done before 4 to 6 months postoperatively.

### PRECAUTIONS

Geistlich Bio-Gide® Shape should only be used by trained dentists or oral surgeons.

In case of membrane exposure during the healing phase the resorption time may be accelerated.

A postoperative maintenance phase can help to ensure long-term therapeutic success.

If endosseous implants are involved, the membrane should be used only in combination with a stable implant and not in lieu of achieving primary implant stability. The long-term safety and effectiveness of the membrane in conjunction with bone-filling materials has not yet been established.

Absolute stability of the membrane is important for guided bone regeneration and is a vital condition for therapeutic success, and the smallest movement on the tissue underneath is to be avoided.

Effect on pregnant or lactating women is not known.

Geistlich Bio-Gide® Shape should be used with special caution in patients with

- Acute or chronic infection (e.g., osteomyelitis) at the surgical site
- Uncontrolled metabolic diseases (e.g., diabetes, osteomalacia, thyroid disorder)
- Prolonged corticosteroid therapy
- Autoimmune diseases
- Radiotherapy
- Heavy smoking

### ADVERSE REACTIONS

As Geistlich Bio-Gide® Shape is a collagen product allergic reactions may not be totally excluded. Possible complications which may occur with any surgery include swelling at the surgical site, flap sloughing, bleeding, dehiscence, hematoma, increased sensitivity and pain, bone loss, redness, and local inflammation.

### STORAGE AND HANDLING

If the sterile packaging is damaged or opened, the product must not be used. The contents of the double blister is designed for single use only.

Do not resterilize Geistlich Bio-Gide® Shape. Resterilization is not possible. Discard any unused material. Do not use after expiry date.

Store at controlled room temperature 15°–25°C (59°–77°F) and in a dry place.

The membrane should be handled using sterile gloves or sterile instruments.

### PRESENTATION

Geistlich Bio-Gide® Shape is packed in sterile double blisters.

| Re-order # | Size   |
|------------|--|
| 500377     | 1 Geistlich Bio-Gide® Shape membrane, 14 x 24 mm |

### Distributed by:

Geistlich Pharma North America Inc.  
Princeton, New Jersey 08540  
1-(855) 799-5500  
[www.geistlich-na.com](http://www.geistlich-na.com)

### Manufacturer:

Geistlich Pharma AG  
CH-6110 Wolhusen  
Switzerland  
Made in Switzerland

CAUTION: Federal law restricts this device to sale by or on the order of a dentist or physician.  
Package insert: [www.geistlich-na.com](http://www.geistlich-na.com)

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